

	QuXAT Quality Management System				Document No:	QUXAT/RF/67
	Title	Discard Autoclave – Validation Register				Page No:

Autoclave Name:			Equipment No:			Sr. No:				
S. No.	Date of Validation	Autoclave Start Time/End Time	Biological Indicator LOT Number	Biological Indicator Expiry Date	Autoclave Temperature & Pressure	Result (Color change) – Yes/No	Biological Indicator (Pass/Fail)	Validation Report Date	Report Reviewed by (Quality Manager)	Validation Approved by (Lab Director)

Discard Autoclave validation process using Biological Indicators is conducted every week.

Refer to Discard Autoclave General SOP & Validation SOP

121°C is a standard temperature for autoclave operation, and generally achieved when chamber pressure is 15-16 psi

Issue No: 1	Issue Date: 1.1.2024	Amend No: 00	Amend Date:	Copy No: 01
Prepared/Issued by:	Quality Manager	Approved by	Laboratory Director	